

DentalMarketplace — Ultrasonic Activator Market Comparison

Eighteenth Ultra X vs. commonly referenced endodontic irrigant activation devices. Competitor prices reflect published U.S. list/distributor pricing as of mid-2026 and may vary by dealer and configuration.

Model	Brand	List Price	Technology	Frequency	Cordless	Weight	Tier
Ultra X	Eighteenth	\$458.99 MAP	Ultrasonic	45 kHz	Yes	95 g	Value
EndoUltra	Vista Apex	\$584.99+	Ultrasonic	40 kHz	Yes	N/A	Premium
EndoActivator	Dentsply Sirona	~\$500–605	Sonic (not ultrasonic)	N/A (sonic)	Yes	N/A	Mid

Why Dentists Choose the Ultra X

Ultra X delivers true ultrasonic activation (45 kHz) — the same underlying technology as the premium EndoUltra, and a meaningfully different mechanism than sonic-only devices like the EndoActivator — at a lower price point than either. For practices that want genuine acoustic streaming and cavitation for irrigant activation, GP removal, and instrument retrieval without paying a premium-brand price, Ultra X is a strong value-tier alternative.

Competitor prices are approximate, sourced from published distributor/dealer list pricing (e.g., PureLife Dental, DentrealStore) as of mid-2026 and are provided for general market positioning only; actual pricing varies by dealer, promotion, and configuration. Specifications sourced from manufacturer documentation and published dealer/distributor listings. Full source list available in the Product Factory QA Report.