

Fast-Pack Pro Heating Tips — Size Buying Guide

Obturation handpiece heating tips are proprietary to each device and are not interchangeable across brands — there is no cross-brand “alternative” tip that fits the Fast-Pack Pro. This guide instead compares the three available sizes to help you choose the right one (or stock a full set).

Spec	Small	Medium	Large
Color	Black	Yellow	Blue
Taper size	40/0.025	50/0.05	60/0.06
Best for	Narrower canals, precision	Versatile mid-size	Wider canals
Price (DentalMarketplace)	\$43.99 MAP	\$43.99 MAP	\$43.99 MAP

Why Stock a Full Set

Canal anatomy varies case to case, and having all three sizes on hand lets clinicians match tip diameter to the canal in front of them rather than compromising on a single option. Keeping a spare of whichever size you use most also means your Fast-Pack Pro stays in rotation instead of waiting on the autoclave.

Important: these tips are compatible with the Fast-Pack Pro only, not the original (non-Pro) Fast-Pack device. Full source list available in the Product Factory QA Report.