



Whicam Story3 CS – Wireless Intraoral Camera

Enjoy Wireless & Auto Focus

Wireless Autofocus Intraoral Camera with Mouse

1 Features of Whicam Story3 CS

-  **High speed digital wireless**
High speed digital wireless provide a clearer live image without interference.
-  **Auto focus**
Real Autofocus function, semi-focus also available.
-  **HD image**
High definition, 1280x720 Image
-  **Gyro mouse built-in**
Useable as a mouse.
-  **Adjustable head angle**
Three adjustable angles to reduce arm fatigue.
-  **280° twist head**
Easy to maneuver in patients mouth.
-  **Larger battery capacity**
User replaceable and longer battery life, adding more hours of usage.
-  **Software compatibility**
Integrates with almost all softwares or we will provide a free imaging software.
-  **OLED screen**
Bright and clear OLED screen indicates various device status.
-  **Adjustable LED brightness**
LED intensity adjustment button, you can adjust brightness to three intensity levels.



2 Wireless Receiver

- USB power operation
It doesn't need other electric power supply.
The receiver works with USB electric power only.
- Compact design
Compact receiver saves space.
- Unlimited receiver channels
Digital wireless type supports unlimited channels.
User can use many receivers with only one camera.

4 Minimum System Requirements

PC system:	2.0 Ghz CPU
OS:	Windows 7 or 8 Windows XP Not HIPAA
USB port:	USB 2.0
RAM:	4 GB RAM

3 Specifications

Camera (included head)	Operation voltage	3.7V
	Measurement	Approx. 266 x 31 x 37 mm
	Weight	Approx. 140g
	Resolution	HD 720p
	Sensor/Lens type	CMOS / Liquid lens
Charger	Measurement	Approx. 169 x 148 x 89 mm
	Weight	Approx. 170 g
	Operation voltage	5V
Disposable cover	Dimensions	Approx. 32 x 100 mm
	Quantity	Approx. 100 pcs by weight x 2ea
	Raw material	Polyethylene
Adapter	Input voltage	AC 100 ~ 240 V
	Output voltage	DC 5V 2A
	UL number	E333054

5 Integration

Integrates with almost ALL software's; or we provide a free Imaging software.

