

Which iSonic Accessory Fits Your Ultrasonic Cleaner?

Model-by-model selection guide | [DentalMarketplace](#)

Stainless Steel Baskets — one per cleaner model

Basket	Fits These Cleaners	Tank Size	Key Difference	Price
SOB4831A	P4831, P4831-NH	9.8" x 6.9" x 3.1" (3 L)	Ships included with the P4831-NH.	\$99.95
SOB4862A	P4861, P4862	12.7" x 8.1" x 3.9" (6 L)	Fits 6" x 8" cassettes.	\$109.95
SOB6.2A	CS6.2	Not published	Standard basket for the CS6.2 tank.	\$109.95
SOB4875A	P4875, P4875II, P4875II-4T, P4875-NH, P4875II-NH, P4875II-4T-NH	13" x 9.4" x 4.1" (7.5 L)	Fits full-size 7 3/4" x 11" cassettes. Six models.	\$119.95
SOB8.0A	CS8.0, P4875-NH-8	13" x 9.4" x 4.1" (7.5 L)	8.3" WIDE basket — the only one that accepts Hu-Friedly Large Thin Ortho Cassettes (8.2" x 11").	\$129.95

Removable Cleaning Chambers — isolate solutions from your tank

Chamber	Fits These Cleaners	Notes	Price
IT4831A	P4831 series	Sized for the P4831 tank.	\$24.95
IT4862A	P4861, P4862, CS6.2, P4890II	Widest compatibility of the chamber range — four models.	\$24.95
IT4875A	P4875, P4875II, CS8.0	Inside approx. 4.75" x 8.75" x 3.25". Two fit side by side in the CS8.0.	\$24.95
IT8.0A	CS8.0, P4875, P4875II	Large single chamber — more working volume than the IT4875A.	\$29.95

Stainless vs. Plastic Baskets

	Stainless Steel Basket	Plastic Basket
Ultrasonic transmission	Absorbs less energy — more cavitation reaches instruments	Absorbs more ultrasonic energy
Cleaning result	Better, per manufacturer	Reduced
Durability	More durable over repeated cycles	Shorter service life
Cleanliness	Easier to keep clean and disinfect	Can retain residue
Noise	Slightly louder in operation	Quieter

Sources: iSonic Inc. manufacturer product pages and collection listings (isonicinc.com) for all fitment, dimensions, pricing, and the stainless-vs-plastic performance rationale. Prices shown are DentalMarketplace list prices and match iSonic published retail. iSonic does not publish standalone basket dimensions; tank sizes are shown so you can confirm the right unit. Baskets and chambers are model-specific and are not interchangeable — confirm your cleaner's model number before ordering. Accessory specifications are subject to manufacturer change.