

Software Name	ensor Support	Single or Progression	Notes
Adstra			
AirTechniques VisionX	Yes		Not to sure about image progression yet
Akitu One	Yes	Progression	TWAIN V1.0.4.0 works grade with Akitu One, Great software
Apteryx DTX Studio	Yes	Progression	Works with all versions of the TWAIN driver, simple set
Apteryx XrayVision	Yes	Progression	versions of this software might possibly have issues V3.12.0.44.
Apteryx XVCapture	Yes	Progression	XVCapture
Cadi Optimum	Yes	Progression	to CADI
CareStream/Softdent	Yes	Single / No Template	
ClearDent	Yes	Progression	driver under "Configuration". However, it usually does crash making it
CliniView Instrumentarium			
ClioSoft/Sotalmage	Yes	Progression	
Cloud9 Software			
CurveDental/Curve Hero	Yes	Progression	templates are already premade. we might just need to change the
DentalEMR			
Dentrix Patient Management	No		DO NOT get confused with Dentrix Ascend
DentrixAscend	Yes	Progression	you locate the tooth in the hidden icons section and right click. Press
DentiMax	Yes	Progression	Templates have to be updated to allow use of Twain.
Dexis	Yes	Progression	However, Dexis support seems to no longer want to help offices with
Dolphin Imaging	Yes	Progression	Plug N play
DoveTail			
Easy Dental Patient Management	No		
EVAsoft	Yes	Progression	
Harmony Jazz	Yes	Progression	Progression requires editing layouts.
iDentalSoft	Yes	Progression	Use D1.0.1.3 driver.
i-Dixel	Yes	Progression	
MediaDent			
MiPACS	Yes	Progression	
MOGO			
OpenDental	Yes	Progression	yet. This is what OpenDental calls progressions.
Oryx Imaging	Yes	Progression	Newest Version of the TWAIN Driver v1.0.4.0 has been running into
Patient Gallery	No		
Patterson Imaging/EagleSoft	Yes	Single	acquire once, take all of the images they need, and then close the
Planmeca Dimaxis Pro			
Planmeca Romexis	Yes	Progression	Planmeca Romexis 6.0.1.812 this version has progression, Twain Driver is in the drop down menu and does capture images, no templates to edit our driver
PortView	No		
Progeny	No		The only Progeny version I've worked with crashed everytime the "Acquire Twain" button was pressed. This happened on every computer.
ProImage			
QuickImageDIS	Yes	Progression	
Sigma Digital			
Sirona Schick	Yes	Single	Using TWAIN drivers with this software is a bit of a hassle for an assistant.  They have to open a template, Ctrl + A to select every template box, right click and select "Acquisition Mode > Scanner", and then click each box every time they want to acquire an image.  I wouldn't recommend the sensor to anyone with this software. I imagine all of those steps needing to be done each time they want to do a progression is going to have them sending the sensor back after some time.
Sirona Sidexis XG	Yes	Single	
SmartDent	Yes	Progression	
SoPro Imaging	Yes	Progression	
Sota Imaging	Yes	Progression	SOTA works great
ProfSuni	No		Recently Prof.Suni
Tab32	Yes	Progression	
TDO	Yes	Single	Few times installed seems to run across some issues but we can get it to work single image, no progression we've contacted TDO to confirm.
Teledentix	Yes	Progression	
TigerView	Yes	Progression	
Tracker	Yes	Progression	Works with porgression, seems like an older software but works fine with newer driver
Triana	Yes	Progression	
UmbieDentalCare			
Vatech EzDent-i	Yes	Single	Newer versions of EzDent-i seem to support single image progression with our driver now
Visix Imaging	Yes	Progression	Visix Imaging only supports one sensor at a time.  Go into the registry and change both instances for the sensor to "TWAIN2 Sensor".  1: Computer\HKEY_USERS\...\VirtualStore\MACHINE\SOFTWARE\WOW6432Node\Air Techniques\Visix\2\Xray\Accent\Options > TwainSourceN... > Change from current sensor to "TWAIN2 Sensor"  2: Computer\HKEY_LOCAL_MACHINE\SOFTWARE\WOW6432Node\Air Techniques\Visix\2\Xray\Accent\Options > TwainSourceN... > Change from current sensor to "TWAIN2 Sensor"  "TWAIN2 Sensor" is going to be whatever name Thwacker/Twain Test gives for our sensor driver.
VixWin Platinum	Yes	Single	Copy twain_32 files to VixWin folder.
XDR	Yes	Progression	XDR works with both versions of the TWAIN Drivers, We would need to install our TWAIN Driver go into the XDR settings and add our device into the TWAIN driver section. our sensors are recognized as TWAIN Drivers not dircetly intragrated into the sensor area. Images are still the same but might have less capabiltity to change and modify images within XDR Settings. In conclusion XDR is one of the better softwares our sensors works with our sensorsunfotuntley they are are recognized under a TWAIN Driver so we might lose some features but if they are okay with that aspect its a great software.